


GUTR Face Shield Protection for the Frontlines



Dedicated to Innovation



Features

- **COMFORT=FOCUS:** GUTR sweatbands are made with a soft, flexible polymer designed to feel comfortable against your skin over long periods of use. They feature an adjustable elastic closure with fine levels of adjustments to ensure a comfortable fit.
- **CLEAR VISION:** GUTR sweatbands capture and redirect sweat from the top of your head and forehead away from your eyes and glasses. Your vision will stay clearer. You won't need to wipe sweat away from your forehead.
- **FACE PROTECTION:** Avoiding flying fluids in your eyes, nose, ears, and mouth is critical for personal safety. The shield is directionally-adjustable for preferred operational positions—and can be raised fully for complete facial access without needing to remove sweatband.
- **AVOID WASTE:** Headbands can be sanitized and reused. Visors can be replaced as needed.
- Can be worn with glasses and mask.
- Made in USA

Specifications

- 0.015" Clear Poly Face Shield - 13" W X 9" Tall
- GUTR sweatbands are made from medical-grade thermoplastic elastomer (TPE)
- Adjustable elastic Strap (1 size fits all)
- Metal fasteners and snaps are nickel-free metal
- Remove Safety Film prior to first use





**PPE BUY
DIRECT**

WWW.PPEBUYDIRECT.COM
sales@ppebuydirect.com



**Providing Low Cost
Sourcing of PPE
At the Best Price!**



Replacing the visor is easy. Just snap-off the old face shield and snap-on a new replacement.

PPE BUY DIRECT is a Division of **INVAC SYSTEMS INC**, a Global Supplier of Equipment, Spare Parts, Service, and Process for the Solar, Semiconductor, MEMS, and Healthcare Industries.



**PPE BUY
DIRECT**

INVAC SYSTEMS



**PPE DIVISION
US SPARE PARTS
TECHNICAL SALES
Dallas - Texas**

**211 N. Dooley Street
Grapevine, TX 76051
(P) +1 817 527 6084**

WWW.PPEBUYDIRECT.COM

