

Level 1 Disposable Surgical Gowns



Dedicated to Innovation

Isolation Gowns Specification

AAMI - Level 1



Product	
Name	Disposable Isolation gown
Material	Medical grade polyethylene 20 grams
Color	Blue
Packaging	1 gown per bag / 100 per case

Specifications	
Performance	Disposable liquid and fluid repellent product Easy tear pattern to be disposed of with ease and hygiene Back adjustment straps Sleeve finished with punch hole for thumb
Measurements	Large: 42.5in x 35.5in x 40in (length x width x sleeve) X-Large: 46.5in x 37.5in x 40in (length x width x sleeve) Length measured at front middle of neckline to bottom Width or circumference measured at waist Sleeve measured from the front center of the neckline
Storage	Store in a fresh and dry place
Purpose	Disposable


MADE IN MEXICO



WWW.PPEBUYDIRECT.COM
sales@ppebuydirect.com



**Providing Low Cost
Sourcing of PPE
At the Best Price!**



PPE BUY DIRECT is a Division of **INVAC SYSTEMS INC**, a Global Supplier of Equipment, Spare Parts, Service, and Process for the Solar, Semiconductor, MEMS, and Healthcare Industries.



**PPE BUY
DIRECT**

INVAC SYSTEMS



**PPE DIVISION
US SPARE PARTS
TECHNICAL SALES
Dallas - Texas**

**211 N. Dooley Street
Grapevine, TX 76051
(P) +1 817 527 6084**

WWW.PPEBUYDIRECT.COM

